


No-D NMR as low-cost and high-throughput measurements

No-D NMR is a technique for measuring NMR signals of organic compounds without deuterated solvents but with normal solvents.

No-D NMR does not require deuterated solvents and so can save costs. Also, the process of dissolving samples in deuterated solvents is omitted, or reaction solution or distillation solution can be measured, as it is, in no-D NMR. Thereby, it can save total experimental time and can be applied to the samples unstable in preparation.


Processes of no-D NMR measurements

1. Resolution adjustments by ^1H -selective FG shimming for normal solvent signals.
 2. Automatic detection and suppression of normal solvent signals.
 3. Automatic referencing using normal solvent signals.
- * NMR lock is not executed.



Automation tool

All the processes automatically run in ECA/ECX series.



Method menu for no-D NMR

Click

1D 2D suite no d


Samples are dissolved in chloroform.

Resolution adjustment, solvent peak suppression, and reference setting are automatically executed.

< Reference >

No-D NMR Spectroscopy: A Simple Yet Powerful Method for Analyzing Reaction and Reagent Solutions


Thomas R.Hoye, Brian M.Eklow, Troy D. Ryba, Mikhail Voloshin, and Letitia J.Yao
Org.Lett., 6.953-956 2004



Results from "1D 2D suite no d"

< Resolution adjustment by ^1H -selective FG shimming >

Resolution can be adjusted by FG shimming for a selected ^1H signal of the sample, such as a normal solvent signal. The result of ^1H -selective FG shimming is shown below.



Before FG shimming


After FG shimming

< Automatic detection and suppression of a normal solvent signal >

An automatic measurement of no-D NMR executes a scout scan to find a solvent signal to be suppressed.



WET solvent suppression is performed by setting a solvent peak position detected automatically (left).



WET can suppress ^{13}C satellites as well as a main peak (right).

- Note -

When solvent peaks are suppressed in no-D NMR, overlapped or adjacent peaks may be lost or reduced. In such cases, solvents should be chosen properly.