


# H2BC & Edited H2BC


Selective observation of  ${}^2J_{CH}$

Long-range correlation peaks in HMBC measurements appear depending on the magnetization transfer time  $1/2 {}^nJ_{CH}$ . However, the value of  ${}^nJ_{CH}$  is not determined by the bond number, and so the correlation peaks give no information on the bond number. In particular, the close values of  ${}^2J_{CH}$  and  ${}^3J_{CH}$  make it difficult to discriminate the corresponding correlations. H2BC (Heteronuclear 2-Bond Correlation) detects only two-bond away correlations, clarifying the assignments. Furthermore, multiplicity edited H2BC may discriminate attached carbon numbers.




Edited H2BC spectrum of strychnine.

${}^1J_{CH}^{max} = 170$  Hz,  ${}^1J_{CH}^{min} = 125$  Hz,  $y\_points = 128$ ,  $T = 22$  msec




H2BC and HSQC spectra of strychnine.



Spectrometer: JNM-ECX400

By linking correlation peaks in H2BC and HSQC spectra,  ${}^{13}C$  connectivities can be revealed. Thus, the combination of H2BC and HSQC gives a method for determining sequences of  ${}^{13}C$  bonding with  ${}^1H$ .



Pulse diagram of H2BC and Edited H2BC measurements.

$$\tau_1 = 1/2\{^1J_{\text{CH}}^{\text{max}} + 0.07(^1J_{\text{CH}}^{\text{max}} - ^1J_{\text{CH}}^{\text{min}})\}^{-1}, \quad \tau_2 = 1/(^1J_{\text{CH}}^{\text{max}} + ^1J_{\text{CH}}^{\text{min}}),$$

$$\tau_3 = 1/2\{^1J_{\text{CH}}^{\text{max}} - 0.07(^1J_{\text{CH}}^{\text{max}} - ^1J_{\text{CH}}^{\text{min}})\}^{-1}, \quad \epsilon = t(\pi^{\text{H}}),$$

$$\delta = \text{PFG delay}, \quad \delta' = \delta + t(\pi^{\text{H}}), \quad 14 < T < 22 \text{ msec}$$

*Reference*

N.T. Nyberg, J.Ø. Duus and O.W. Sørensen, *J. Am. Chem. Soc.* 127, 6154 (2005)  
 N.T. Nyberg, J.Ø. Duus and O.W. Sørensen, *Magn. Reson. Chem.* 43, 971 (2005)