


## $^{25}\text{Mg}$ solid-state NMR in an ultra-high magnetic field: Structures of bacteriochlorophyll *c* complex, chlorosome

Mg atoms play an important role in biological materials, but  $^{25}\text{Mg}$  NMR measurements are difficult since  $^{25}\text{Mg}$  is low- $\gamma$  nucleus and exhibits quadrupolar broadening in solids. An ultra-high magnetic field may, however, dramatically enhance sensitivity.


This Note introduces a study on stacking structures of bacteriochlorophyll *c* in chlorosome clarified using  $^{25}\text{Mg}$  solid-state NMR at 21.8 T.

\* This is collaborative research with Prof. Y. Koyama of Kansai Gakuin Univ., Prof. H. Nagae of Kobe City Univ. of Foreign Studies, and Prof. T. Shimizu of National Institute of Materials Science.




Chlorosome is a light-harvesting organ of green photosynthesis bacteria, and bacteriochlorophyll *c* (BChl *c*) aggregate in it. Stacking structures of BChl *c* have attracted interest for a long time.


Using a 21.8 T ultra-high field,  $^{25}\text{Mg}$  MAS NMR spectrum of chlorosome was obtained within 80 hours, unprecedentedly short time.



$^{25}\text{Mg}$  spectrum of chlorosome suggests two Mg sites: Large quadrupolar or low symmetry Mg (Site 1) and small quadrupolar or high symmetry Mg (Site 2). Fig. 3 shows a plausible model: Mg-coordinated OH deviates Mg atom of Site 1 from the ring plane, whereas Mg atom of Site 2 stays within the plane because of a ligand water molecule.



Thus, BChl *c* molecules in chlorosome prove to have stacking structures depicted in Fig. 4.



This study demonstrates that ultra-high field NMR system is useful for low- $\gamma$  nuclei in real samples, such as biological materials.

Reference: Y. Kakitani, Y. Koyama, Y. Shimoikeda, T. Nakai, H. Utsumi, T. Shimizu, and H. Nagae, *Biochemistry* **2009**, *48*, 74-86.