


## STMAS probe

*to analyze low-sensitivity half-integer spin nuclei in solid-state high-resolution NMR spectroscopy*

JEOL developed a special probe for realizing STMAS (Satellite Transition Magic Angle Spinning) measurements which give high-resolution NMR spectra for half-integer spin nuclei in solids. STMAS may give much higher spectral sensitivity in comparison with MQMAS (Multiple Quantum Magic Angle Spinning).



### Comparison of sensitivity in STMAS and MQMAS

Spectrometer: JNM-ECA700, Probe: 3.2mm STMAS


Sample: Kaolin ( $Cq = 3.5 \pm 0.4$  MHz) and glass ( $Cq = 6.8 \pm 0.6$  MHz)

Spectral sensitivity is compared for various RF field strengths and pulse sequences: Z-filter 3QMAS (+), 3QMAS-SPAM (■), DQ-STMAS (△), DQF-STMAS (○), and DQF-STMAS-SPAM (●).

Reference: T. Takahashi, K. Kanehashi, Y. Shimoikeda, T. Nemoto, and K. Saito, J. Magn. Reson., 198, 228 (2009).



## Five distinct features of STMAS probe



JEOL STMAS has special mechanism for **accurate magic angle adjustment** within 0.001 degree and **its long-term retention**, required in STMAS measurements. STMAS probe, by combining with MAS Controller, allows **measurements for low-sensitivity samples** and **great reduction of measurement time**.

JEOL STMAS probe was developed in collaboration with Dr. K. Kanehashi and Dr. T. Takahashi of Nippon Steel Corporation.