


## Homonuclear covariance transformation in solid-state NMR!

Covariance (Cov) transformation helps us to reduce the data collection time in homonuclear correlation solid-state NMR experiments such as NOESY and INADEQUATE. In these experiments, large number of  $t_1$  points are usually required to achieve sufficient spectral resolution in the indirect dimension. For example, fully collected Fourier transformed (FT) spectrum of L-histidine monochloride was obtained with 512  $t_1$  points (A). When the number of  $t_1$  points was reduced to 40 (B), significant wiggle was observed, resulting in very low resolution. Cov spectra (C) of 40  $t_1$  points (C) produces a spectrum almost identical to (A). Cov achieves more than 10 times reduction of experimental time!


$^{13}\text{C}$ - $^{13}\text{C}$  PDSD  
correlation  
spectra

ECA500  
3.2mm HXMAS  
20 kHz MAS  
 $^{13}\text{C}_6$ ,  $^{15}\text{N}_3$  L-histidine HCl




(A) FT: 512  $t_1$  points

1/10 time  
without losing  
anything!



(B) FT: 40  $t_1$  points

Covariance!



(C) Cov: 40  $t_1$  points

J. Magn. Reson. 171 (2004) 277-283.

Copyright © 2012 JEOL RESONANCE Inc.

<http://www.jeol.co.jp>  
<http://j-resonance.com>



**JEOL RESONANCE**