


Space-saving carousel-type 24-slot auto sample changer

ASC24

< Features >

- Permit auto change of 24 sample tubes.
- Space-saving installation at the top of a magnet.
- Removable carousel-type container.
- New controller unit attached.



A sample container is removable when loading samples.



Controller unit

< Specifications >

- Main unit with a sample container
 - Number of sample tubes: 24
 - Sample tube requirements:
 - 5mm ϕ or 10mm ϕ tubes for solution NMR with the length of 178 to 203mm
 - (* Please inquire about 3mm ϕ tubes)
 - Dimensions: 325mm x 400mm x 350mm
 - Weight: 12.4kg
 - (mechanical unit 8.2kg, sample container 2.8kg, sample tubes and holders 1.4kg)
- Controller unit
 - Power source: single-phase 100-125V 50/60Hz 1A
 - Dimensions: 280mm x 355mm x 158mm
 - Weight: 4kg
- Required system: JNM-ECA/ECX series with Delta 4.3.6 or later